IPS Empress® Direct Color

[en] Instructions for Use

Stains

Rx ONLY

CE 0123

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

Date information prepared:
2024-06-26 / Rev. 0
web

ivoclar

English

1 Intended use

Intended purpose

Esthetic characterization of direct and indirect anterior and posterior restorations and masking of discolorations.

Patient target group

Patients with permanent teeth.

Intended users

Dentists

Special training

No special training required.

Use

For dental use only.

Description

IPS Empress Direct Color is a light-curing characterization material, which is used to characterize esthetically demanding composite restorations in the direct restorative treatment of anterior and posterior teeth.

IPS Empress Direct Color cures with light in the wavelength range of 400–500 nm.

Indications

- None

Types of restorations:

- Individualized characterization of direct and lab-fabricated composite resin restorations
- Masking of discolored tooth structure
- Characterization of chairside CAD/CAM veneers.

Contraindications

The use of these products is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- If a dry working field cannot be established.
- If the stipulated working technique cannot be applied.
- When used as a restorative without a capping layer.
- Apply the product at ambient temperature. Cool temperatures render the material difficult to extrude.
- The use of different application cannulas can render the material difficult to extrude.
- The syringe may be used for multiple applications, but should be protected by a single-use hygienic sleeve. Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
- ☒ The application cannula of the syringe must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).

Side effects

In rare cases, components of IPS Empress Direct Color may lead to sensitization. The products must not be used in such cases.

Interactions

- Phenolic substances such as eugenol/clove oil inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with IPS Empress Direct Color must be avoided.
- Discoloration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Clinical benefit

- Restoration of esthetics.

Composition

- **Brown shade:**
Bis-GMA, D3MA, UDMA, silicone oxide, TEGDMA, camphorquinone, aminobenzoate
- **Remaining shades:**
Bis-GMA, UDMA, TEGDMA, silicone oxide

Total content of inorganic fillers: 5.0–6.5 vol%.

Particle size of inorganic fillers: between 0.03 µm and 16.3 µm.

2 Application

The product is used in combination with composite resin materials as well as for a special technique to customize veneers fabricated using chairside CAD/CAM. The corresponding pre-treatment with other products should be carried out as specified in the instructions for use of that particular material. The product (IPS Empress Direct Color) can be dispensed from the syringe and applied with a brush or a probe. However, the product can also be applied directly to the restorative material with the thin, pre-bent metal cannulas that are included in the delivery. In order to create very small spots of colour, it is advisable to wipe the cannula with a cellulose pad from time to time to remove any excess material.

Notice! The maximum increment is 0.5 mm for all shades except for brown. The maximum increment for brown is 0.1 mm.

Shade determination

The product is available in the following shades:

Shade	Possible applications
White	<ul style="list-style-type: none">– Imitating enamel stains– Simulating halo effects– Accentuating cusps and ridges
Honey Ochre	<ul style="list-style-type: none">– Emphasizing the dentin core– Highlighting fissures– Imitating individual enamel characteristics
Blue Grey	<ul style="list-style-type: none">– Mimicking translucent enamel areas– Imitating opalescent areas
Brown	<ul style="list-style-type: none">– Imitating discolored fissures, pits and enamel cracks

Isolation

Adequate relative or absolute isolation is required.

A) Use for the characterization of direct restorations

1. Use the composite resin material according to the manufacturer's instructions. When restoring a large surface, IPS Empress Direct Color can be placed during the build-up of the restoration and immediately covered with a layer of composite. Fissures and enamel cracks can be simulated as follows:
 - a) In the appropriate place, carve a notch in the last ("enamel") composite increment.
 - b) Fill the stain material into the notch using the cannula or a suitable instrument.
 - c) Then, cover the notch with the composite resin material, starting from the sides. The stain material rests within the last composite layer and is therefore protected from wear.

Alternative delayed application of the product:

1. Grind in a "characterization pattern" into the restorative material.
2. Apply a 37% phosphoric acid gel and allow the gel to react for 15 s.
3. Rinse off the phosphoric acid gel with water and dry the surface with moisture- and oil-free air.

4. In order to establish a sound bond with the polymerized restorative material, wet the cleaned surface with a light-curing bonding agent.
- Notice!** The maximum increment is 0.5 mm for all shades except for brown. The maximum increment for brown is 0.1 mm.
5. Apply the product with a brush or another suitable instrument.
6. Light-cure the product according to the indicated exposure time and light intensity in Table 1.

 **NOTICE! The instructions for use of the curing light used must be observed.**

7. Layer the surface with the composite resin material and light-cure the material according to the manufacturer's instructions.

B) Use for the characterization of lab-fabricated composite restorations

When layering inlays/onlays, the product can be applied in a manner similar to the process used for direct restorations (section A).

C) Use for the characterization of chairside CAD/CAM restorations

Proceed as follows:

1. Try-in and subsequently adjust the veneers.
2. Clean and condition the inner surfaces of the ceramic veneers according to the manufacturer's instructions.
3. Apply a thin layer of a light-curing bonding agent to the conditioned inner surfaces of the veneers and disperse this adhesive layer to a very thin film with blown air.
- Notice!** Excessively thick shade layers could compromise the fit and the esthetic appearance of the restoration.
- Notice!** Avoid adjustments with rotating instruments as well as additional try-ins after characterization. Otherwise, the bond to the luting material can no longer be ensured.
4. Apply very small amounts of the product in thin layers in order to meet individual characterization requirements.
5. Light-cure the light-curing bonding agent and the stains in accordance to the information in Table 1.
-  **NOTICE! The instructions for use of the curing light must be observed.**
6. Place the veneer according to the instructions of the manufacturer of the luting material.

Application notes

- If necessary, the products can also be applied directly to already polymerized material. If the restoration is already polished, it must be roughened and wetted first with an adhesive before the new product can be applied. (IPS Empress Direct Opal shade and IPS Empress Effect are not suitable for the reconstruction of occlusal surfaces).
- Only cannulas intended for the product must be used. In case other cannulas are used, there is a risk that the cannula may suddenly become loose during application.
- Do not expose the product to intensive light while applying it, as this might lead to premature polymerization.
- To ensure long-lasting shade effects, the product should not be applied as final layer on the surface of the restoration, but, if possible, it should be placed as intermediate layer or in carved notches.
- Do not disinfect syringes with oxidizing disinfection agents.
- Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
- Avoid any contact of the reusable syringe not covered with a protective sleeve with the mouth of the patient.
- Dispose of contaminated syringes.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website www.ivoclar.com
- Explanation of the symbols: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL015JM

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Avoid contact of the uncured product with the skin/mucous membrane or eyes. The unpolymerized product may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity
- Inclusion of air bubbles during the placement of the filling
- Chipping, fractures
- Ingestion of material

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringes immediately after usage. Exposure to light leads to premature polymerization.
- Do not use the product after the indicated date of expiration.
- Expiry date: See information on syringes and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

Table 1

Light intensity	Exposure time
500–900 mW/cm ²	20 s
1,000–1,300 mW/cm ²	10 s

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