

# Product information

## FD 366 sensitive

### disinfection of sensitive surfaces



#### FD 366 sensitive at a glance

- Aldehyde-free, alcohol-based ready-to-use solution for quick-acting disinfection of sensitive surfaces on invasive and non-invasive medical devices in private practices, laboratories and hospitals.
- Above-average material compatibility. Particularly suitable for sensitive surfaces made of plastic/acrylic glass or synthetic leather on treatment units.
- In accordance with the VAH, the exposure time is 1 min. (bactericidal / yeasticidal).
- Broad spectrum of activity: bactericidal, tuberculocidal, yeasticidal, virucidal activity against enveloped viruses (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses.
- Tested in accordance with current VAH methods and European standards.
- Tested in 4-field test in accordance with VAH and EN 16615.
- VAH list, IHO list of disinfectants.
- Especially recommended for use in combination with FD multi wipes/FD multi wipes compact.
- Dries quickly.

#### Properties

FD 366 sensitive from Dürr System Hygiene is an aldehyde-free, alcohol-based ready-to-use solution for quick-acting disinfection of sensitive surfaces on invasive and non-invasive medical devices and medical furnishings in private practices, laboratories and hospitals. FD 366 sensitive is distinguished by fast efficiency and good material compatibility. Disinfected surfaces dry rapidly. FD 366 sensitive is especially recommended for use with FD multi wipes/FD multi wipes compact.

#### Product composition

FD 366 sensitive contains an active substance combination of alcohols and quaternary ammonium cations. 100 g FD 366 sensitive contains 17 g 1-propanol, 0.125 g alkyl-benzyl-dimethylammonium chloride, excipients and water.

#### Antimicrobial effect

FD 366 sensitive has bactericidal, tuberculocidal, yeasticidal, virucidal activity against enveloped viruses (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses. FD 366 sensitive is included in the VAH list and IHO list of disinfectants. Tested under low and high load in accordance with EN 13727, EN 13624, EN 14348, EN 14476, EN 16615, as well as DVV/RKI guidelines.

#### Directions for use

FD 366 sensitive is applied undiluted. Apply FD 366 sensitive to a cloth and fully wet the surfaces. Wipe the surfaces and objects to be disinfected and allow the product to take effect. Disinfect visibly contaminated surfaces with FD 366 sensitive and remove contamination. Disinfect cleaned surfaces again. Make sure that wetting is complete throughout the entire action time. In accordance with VAH, the action time is 1 minute (bactericidal / yeasticidal), 2 minutes for Tb bacteria. Inactivation of viruses takes place within 30 seconds to 5 minutes. Use FD multi wipes/FD multi wipes compact saturated with FD°366 in the same way.

#### Environmental impact

The packaging is made of polyethylene and polypropylene and can thus be recycled or incinerated. For recycling, flush bottles with water.

For further information see safety data sheet.

#### Shelf life

Do not use the product after the expiry date printed on the packaging. When using with FD multi wipes/FD multi wipes compact, 3 months.

# Product information

## FD 366 sensitive

### disinfection of sensitive surfaces

#### Container sizes

10-L canister  
2.5-L bottle  
1-L bottle  
750-mL bottle for Hygowipe Plus

#### Storage

Store in a well-ventilated place. Keep container tightly closed.

#### Accessories

Spray bottle (600 mL), spray nozzle for 1-L bottle, canister spigot, FD multi wipes/FD multi wipes compact

#### General notes

Disinfect the exterior surfaces of medical devices, such as dental curing lights, nozzles, handpieces, and contra-angle handpieces, prior to processing as required in accordance with the manufacturers' instructions. In the case of sensitive surfaces, test for compatibility in an inconspicuous place. For information on fire and explosion protection, see DGV 107-002 (BGR 206). Do not mix FD 366 sensitive with other products. Changes in product colour/odour may occur, especially when stored in sunlight. However, these

changes in colour and/or odour have no impact on the disinfection efficacy of the product. Use exclusively by qualified personnel. All serious events that occur in association with the product must be immediately reported to the manufacturer and the competent authority of your member state.

#### Hazard warnings

FD 366 sensitive is classified and labelled according to the CLP Regulation. For further information see safety data sheet.  
Flash point: 33°C

#### Independent expert opinions – in-house investigations

The expert opinions are available upon request.



Application		Contact time
Recommendation for use for surface disinfection in accordance with VAH	Bactericidal / yeasticidal	1 min.

Type of activity	Tested according to	Contact time
<b>Bacteria and yeasts</b>		
Bactericidal / yeasticidal	EN 13727, EN 13624, EN 16615	1 min.
Tuberculocidal	EN 14348	2 min.
<b>Viruses</b>		
Virucidal activity against enveloped viruses (enveloped viruses including HBV, HCV, HIV as well as coronaviruses)	EN 14476	30 sec.
	DVV / RKI guideline	30 sec.
Noroviruses	EN 14476	5 min. <sup>1)</sup>

1) Testing under low burden.

#### Sales

Dürr Dental SE  
Höpfigheimer Str. 17  
74321 Bietigheim-Bissingen  
Germany  
www.duerddental.com  
info@duerrdental.com

#### Manufacturer

orochemie GmbH + Co. KG  
Max-Planck-Straße 27  
70806 Kornwestheim  
Germany  
www.orochemie.de  
info@orochemie.com

